 **BEDIENUNGSANLEITUNG**
BODY SHAPE

 **MANUAL**
BODY SHAPE

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CONGRATULATIONS!

With the purchase of the **BODY SHAPE** you have proved consciousness for your health.

To use the advantages of this massage tool for a long time, please read the operation manual thoroughly for the correct methods of usage.

We wish you much joy with your personal **BODY SHAPE**.

SAFETY INSTRUCTIONS

Please read the safety instructions very carefully before using this massage tool. Should you have some of the following health states, please take advice by the current doctor before you take this massage chair in usage.

If you:

- are suffering acute illness
 - suffer a heart illness
 - suffer Osteoporoses
 - are pregnant
 - have a malicious tumour
 - difficulties in perception
 - are getting medical treatment
 - are physically handicapped
-
- If aches arise during the use of the **BODY SHAPE** visit a doctor please immediately.
 - In order to avoid an electric shock or injuries, do not permit your children use the **BODY SHAPE** unsupervised.
 - Do not use the **BODY SHAPE** at places where water and maximum temperatures are exposed.
 - Do not use damaged plugs, cords or loose versions. If plugs or cords are damaged these must be replaced by the manufacturer, qualified personnel or by the service representative.
 - Avoid pulling and rotating the **BODY SHAPE**. In addition, avoid placing heavy articles on it.

- Making sure after long use that the connection is separated in order to avoid the danger of an electric shock, short-circuit or fire.
- With non function of the **BODY SHAPE** please disconnect immediately.
- Avoid touching the plug with wet hands.
- Do not use the **BODY SHAPE** directly after the meal; it could cause nausea.
- The massage time should not exceed the recommended time of 30 minutes, in order to avoid an excessive stimulation of the musculature and the nerves.

MEDICAL REFERENCES

If your physician/therapist did not expressly forbid the back massage, you can enjoy the massage without doubts.

The massage strength increases if the **BODY SHAPE** is brought into a lying position. It is recommended to start massage in a seat position.

The massage should never be felt as unpleasant. Exaggerated massage can also damage the body.

Pregnant women, patients and children should not use the **BODY SHAPE**. With injuries of the spinal column or after operations a massage can be dangerous. In case of doubt, please speak with your doctor.

CARE

You have acquired a high-quality product, which requires the appropriate care. Avoid direct sun and high humidity.

TERMS OF WARRANTY

The dealer grants a legally prescribed guarantee time on the product. The warranty starts on the day of purchase. The purchase date is to be proven with the purchase invoice.

The dealer will repair defect products with determined errors free of charge and/or will replace it in Germany. Wearing parts e.g. material covers are not included. A condition for this is the entry of the defective product with purchase invoice before expiration of the guarantee period.

The warranty claim will be cancelled if a statement of an error, e.g. by outside effect or as consequence of repair or change which was not made by the manufacturer or an authorized contract dealer appeared.


The warranty granted by the dealer is limited to repair and/or exchange of the product. In the context of this warranty the manufacturer and/or dealer, transfers no adhesion and is not responsible for damage which results from ignoring the operating instructions and/or abuse of the product.

If it turns out with examination of the product by the dealer that the warranty demand concerns an error not included by the warranty or that the warranty term ran out the costs of examination and repair is to be carried by the customer.

PRODUCT FEATURES

- 10 vibration motors
- pleasant heat function
- 5 massage modes
- adjustable intensity
- 4 massage zones, each zone can be used separately

BY THE REMOTE CONTROL YOU CAN PILOT THE FOLLOWING FUNCTIONS:

TAKE IN USE

Examine before start-up whether the mains voltage agrees with the tension indicated on the power supply plug.

If not already done connect the mains cable with the back part of the massage chair and plug the power supply into the plug socket.

To start-up the **BODY SHAPE** press ON/OFF.

selection of the massage mode

Choose between 5 massage modes. Press the key „MODE F1 -F5“ to make a choice. The LED-lights show which massage mode is activated.

selection of the massage zone

Choose between 4 massage zones. Each zone can be used separately: neck, upper back, lumbar and legs. Press the adequate key to make a choice.

selection of the massage intensity

Choose between 3 massage intensity modes: low, middle and high.

Press the key „Intensity“ to make your choice. The LED-lights L/M/H show the chosen intensity mode.

heat function

The **BODY SHAPE** has also a heat function. Press the key „Heat“ to activate this function.

preprogrammed timer function

The **BODY SHAPE** has a preprogrammed timer function. This function switches off the **BODY SHAPE** after:

- 15 minutes, when the heat function is activated.
- 30 minutes without heat function.